

Aero Force - Metal Soccer Guard System : Sizing Chart

Effective as of 8/31/2007 (Subject to Change)

Guard Width	PowerChair Frame Width (inches)					Guard Width Type
	14"	15"	16"	17"	18"	
PowerChair Make & Model	Small	Small	Small / Std	Std	Std	
	12	6	0 / 12	6	0	# of 1/4" Spacers
Invacare Arrow (Rear-Wheel Drive)	Small	Small	Small / Std	Std	Std	Guard Width Type
	12	6	0 / 12	6	0	# of 1/4" Spacers
Invacare Ranger (Rear-Wheel Drive)	Small	Small	Small / Std	Std	Std	Guard Width Type
	12	6	0 / 12	6	0	# of 1/4" Spacers
Invacare Torque (Rear-Wheel Drive)	Small	Small	Small / Std	Std	Std	Guard Width Type
	12	6	0 / 12	6	0	# of 1/4" Spacers
Quickie P-200	Small	Small	Small / Std	Std	Std	Guard Width Type
	12	6	0 / 12	6	0	# of 1/4" Spacers
Quickie P-220/222	Small	Small	Small / Std	Std	Std	Guard Width Type
	12	6	0 / 12	6	0	# of 1/4" Spacers
Quickie P-320	Small	Small	Small / Std	Std	Std	Guard Width Type
	12	6	0 / 12	6	0	# of 1/4" Spacers
Quickie S-626/646	Small	Small	Small / Std	Std	Std	Guard Width Type
	12	6	0 / 12	6	0	# of 1/4" Spacers
Boulder Plus H-Frame	Small	Small	Small / Std	Std	Std	Guard Width Type
	12	6	0 / 12	6	0	# of 1/4" Spacers

Small = Small Guard Width (Internal Dimension Width is ~17")

Std = Standard Guard Width (Internal Dimension Width is ~19")

Guard Length	Distance from Reversed Front Tires to Toes	
	<= 12"	13" - 15"
Guard Length Type	Std	Long
Guard Length	12"	16"
Total Length from Reversed Front Tires to Edge of Guard	13"	17"

Std = Standard Guard Length (Outside Dimension Length is ~12")

Long = Long Guard Length (Outside Dimension Length is ~16")

=====